


Arginase 2 (Arg2) Rabbit Pab For Wb, Elisa - P78540

Item Number: CM0029867



Introduction

Arginase 2 (ARG2) rabbit polyclonal antibody for WB and ELISA; detects human, mouse, rat. Immunogen corresponds to amino acids 270-354. Useful for arginine metabolism, immune regulation, and mitochondrial research.

[Learn More](#)

Parameter	Value
Product Name	Arginase 2 (ARG2) Rabbit pAb
Remarks/Alias	ARG2; Arginase 2 (ARG2)
Species	Human
GeneID (human)	384
GeneID	384
Immunogen	Recombinant protein corresponding to amino acids 270-354 of human Arginase 2 (ARG2) (NP_001163.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P78540
Protein Weight	39kDa
Shipping	Ice bag