

Arhgap30 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q7Z6I6

Item Number: CM0022792

Introduction

ARHGAP30 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA in human and mouse. Targets Rho GTPase-activating protein 30, a GAP for RAC1 and RHOA. Recombinant immunogen, rabbit source.

[Learn More](#)

Parameter	Value
Product Name	ARHGAP30 Rabbit pAb
Remarks/Alias	ARHGAP30
Species	Human
Gene ID (Human)	257106
Gene ID	257106
Immunogen	Recombinant protein corresponding to amino acids 485-585 of human ARHGAP30 (NP_859071.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q7Z6I6
Protein Weight	119kDa
Shipping	Ice bag