

Arhgdib Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P52566

Item Number: CM0026016



Introduction

High-quality rabbit polyclonal antibody against human ARHGDIB (RhoGDI2/Ly-GDI). Validated for WB, IHC-P, and ELISA. Detects target in human and mouse samples. Ideal for studying Rho GTPase signaling and actin cytoskeleton dynamics.

[Learn More](#)

Parameter	Value
Product Name	ARHGDIB Rabbit pAb
Remarks/Alias	D4; GDIA2; GDID4; LYGDI; Ly-GDI; RAP1GN1; RhoGDI2; ARHGDIB
Species	Human
GeneID (Human)	397
GeneID	397
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-201 of human ARHGDIB (NP_001166.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P52566
Protein Weight	23kDa
Shipping	Ice bag