

Arms2 Rabbit Polyclonal Antibody For Wb And Elisa - P0C7Q2

Item Number: CM0018136

Introduction

High-quality ARMS2 Rabbit Polyclonal Antibody validated for WB and ELISA, suitable for investigating age-related macular degeneration and ARMS2 protein expression in human retina and placenta. Cross-reacts with mouse.

[Learn More](#)

Parameter	Value
Product Name	ARMS2 Rabbit pAb
Remarks/Alias	ARMD8; ARMS2
Species	Human
Gene ID (Human)	387715
Gene ID	387715
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-107 of human ARMS2 (NP_001093137.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P0C7Q2
Protein Weight	11kDa
Shipping	Ice bag