

Arr3 Rabbit Pab For Wb, Elisa - P36575

Item Number: CM0015598

Introduction

Rabbit polyclonal antibody against human ARR3 (cone arrestin). Validated for Western blot and ELISA, cross-reacts with mouse and rat. Used in retinal signal transduction research.

[Learn More](#)

Parameter	Value
Product Name	ARR3 Rabbit pAb
Remarks/Alias	ARRX; cArr; MYP26; ARR3
Species	Human
GeneID	407
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human ARR3 (NP_004303.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P36575
Protein Weight	43kDa
Shipping	Ice bag