

Asb17 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Wxj9

Item Number: CM0024179

Introduction

Anti-ASB17 Rabbit Polyclonal Antibody validated for WB and ELISA applications, with cross reactivity to Human, Mouse, and Rat. Suitable for research on ubiquitin-proteasome pathway in testis.

[Learn More](#)

Parameter	Value
Product Name	ASB17 Rabbit pAb
Remarks/Alias	Asb-17; ASB17
Species	Human
GeneID (Human)	127247
GeneID	127247
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-130 of human ASB17 (NP_543144.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WXJ9
Protein Weight	34kDa
Shipping	Ice bag