

Ascl1 Rabbit Polyclonal Antibody For Wb And Elisa - P50553

Item Number: CM0017556

Introduction

Rabbit polyclonal antibody targeting human ASCL1 (Achaete-scute homolog 1), a transcription factor crucial for neuronal differentiation. Validated for Western blot (WB) and ELISA; cross-reacts with mouse and rat. Suitable for neurodevelopmental research.

[Learn More](#)

Parameter	Value
Product Name	Ascl1 Rabbit pAb
Remarks/Alias	ASH1; HASH1; MASH1; bHLHa46; Ascl1
Species	Human
GeneID (human)	429
GeneID	429
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 150-236 of human Ascl1 (NP_004307.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P50553
Protein Weight	25kDa
Shipping	Ice bag