


Asc/Tms1 Rabbit Monoclonal Antibody For Wb, If, Ihc, Elisa - Q9Epb4

Item Number: CM0006294



Introduction

Recombinant rabbit monoclonal antibody against mouse ASC/TMS1 (Pycard), validated for WB, IF, IHC, and ELISA. Detects human, mouse, and rat ASC; ideal for inflammasome and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	ASC/TMS1 Rabbit mAb
Remarks/Alias	Asc; TNS1; masc; CARD5; TMS-1; 9130417A21Rik; ASC/TMS1
Species	Mouse
GeneID (Human)	29108
GeneID	66824
Immunogen	Recombinant protein corresponding to amino acids 1-93 of mouse ASC/TMS1 (NP_075747.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9EPB4
Protein Weight	21kDa
Shipping	Ice bag