


Asc/Tms1 Rabbit Polyclonal Antibody For Wb, If, Ihc-P, Elisa - Q9ULZ3

Item Number: CM0015998



Introduction

Rabbit polyclonal antibody against ASC/TMS1 (PYCARD). Confirmed for Western blot, immunofluorescence, IHC-P, and ELISA. Detects human, mouse, and rat ASC/TMS1. For apoptosis, inflammasome, and innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	ASC/TMS1 Rabbit pAb
Remarks/Alias	ASC; TMS; TMS1; CARD5; TMS-1; ASC/TMS1
Species	Human
Gene ID (Human)	29108
Gene ID	29108
Immunogen	Synthetic peptide corresponding to amino acids 50-195 of human ASC/TMS1 (NP_037390.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9ULZ3
Protein Weight	22kDa
Shipping	Ice bag