

Asmt Rabbit Polyclonal Antibody For Wb, Elisa - P46597

Item Number: CM0018238



Introduction

Rabbit polyclonal antibody targeting human ASMT (acetylserotonin O-methyltransferase), for WB and ELISA applications. Ideal for investigating melatonin biosynthesis pathways. Cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	ASMT Rabbit pAb
Remarks/Alias	ASMTY, HIOMT, HIOMTY, ASMT
Species	Human
GeneID (Human)	438
GeneID	438
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 76-298 of human ASMT (NP_001164510.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P46597
Protein Weight	38kDa
Shipping	Ice bag