

Aspm Rabbit Pab - Gene Id: 259266

Item Number: CM0022595

Introduction

Rabbit polyclonal antibody targeting ASPM (abnormal spindle-like microcephaly-associated protein). Validated for WB, IF/ICC, ELISA in human, mouse, and rat samples. Suitable for neurodevelopmental and cancer research.

[Learn More](#)

Parameter	Specification
Product Name	ASPM Rabbit pAb
Remarks/Alias	ASP; MCPH5; Calmbp1; ASPM
Gene ID	259266
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 450-550 of human ASPM (NP_060606.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
Shipping	Ice bag