


Ass1 Rabbit Pab - P00966

Item Number: CM0022279



Introduction

Rabbit polyclonal antibody against human ASS1 (argininosuccinate synthase), validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat; suitable for urea cycle and arginine biosynthesis research.

[Learn More](#)

Parameter	Value
Product Name	ASS1 Rabbit pAb
Remarks/Alias	ASS; CTLN1; ASS1
Species	Human
Gene ID (Human)	445
Gene ID	445
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-412 of human ASS1 (NP_446464.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00966
Protein Weight	47kDa
Shipping	Ice bag