


Asymmetric Dimethyl-Histone H3 (R17) Rabbit Polyclonal Antibody - Q16695/P68431 | Wb, Ihc-P, If/Icc, Elisa

Item Number: CM0000384



Introduction

Rabbit polyclonal antibody specific for Histone H3 asymmetrically dimethylated at arginine 17 (R17). Validated for WB, IHC-P, IF/ICC, and ELISA; cross-reacts with human, mouse, and rat. Ideal for epigenetic and chromatin research studies.

[Learn More](#)

Parameter	Value
Product Name	Asymmetric DiMethyl-Histone H3-R17 Rabbit pAb
Remarks/Alias	H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; Asymmetric DiMethyl-Histone H3-R17
Species	Human
GeneID (Human)	H3
GeneID	8290/8350
Immunogen	Synthetic peptide: A synthetic asymmetric dimethylated peptide around R17 of human Histone H3 (NP_003520.1)
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15kDa
Shipping	Ice bag