

Asymmetric Dimethyl-Histone H3-R2 (Q16695/P68431) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa

Item Number: CM0000271



Introduction

Rabbit polyclonal antibody targeting asymmetric dimethylation at arginine 2 of histone H3 (Q16695/P68431). Validated for WB, IF/ICC, ELISA; cross-reacts with human, mouse, rat and other species.

[Learn More](#)

Parameter	Details
Product Name	Asymmetric DiMethyl-Histone H3-R2 Rabbit pAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; Asymmetric DiMethyl-Histone H3-R2
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	A synthetic asymmetric dimethylated peptide around R2 of human histone H3 (NP_003520.1)
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag