

Asymmetric Dimethyl-Histone H3-R2 Rabbit Polyclonal Antibody (Q16695/P68431) For Wb And Elisa

Item Number: CM0000243



Introduction

This rabbit polyclonal antibody specifically recognizes asymmetric dimethylation of Histone H3 at arginine 2 (H3R2me2a). Suitable for western blot and ELISA, it cross-reacts with human, mouse, rat, and other species, supporting epigenetic and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	Asymmetric DiMethyl-Histone H3-R2 Rabbit pAb
Remarks/Alias	H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; Asymmetric DiMethyl-Histone H3-R2
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	A synthetic asymmetric dimethylated peptide around R2 of human Histone H3 (NP_003520.1)
Source	Rabbit
Category	Methylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15kDa
Shipping	Ice bag