


Asymmetric Dimethyl-Histone H3-R26 (Q16695 / P68431) Rabbit Pab For Wb, If/Icc, Elisa

Item Number: CM0000305



Introduction

Rabbit polyclonal antibody targeting asymmetric dimethyl-histone H3 at R26 (Q16695/P68431), suitable for WB, IF/ICC, and ELISA. Validated in Human, Mouse, Rat, and other wide-range predicted species. Ideal for chromatin and epigenetic studies.

[Learn More](#)

Parameter	Value
Product Name	Asymmetric DiMethyl-Histone H3-R26 Rabbit pAb
Remarks/Alias	H3/j; H3C1; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3Fj; H3C10; H3C11; HIST1H3J; Asymmetric DiMethyl-Histone H3-R26
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	Synthetic peptide
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15kDa
Shipping	Ice bag