

Asymmetric Dimethyl-Histone H3-R8 Rabbit Mab (Q16695/P68431) For Epigenetics Research

Item Number: CM0000249



Introduction

Rabbit monoclonal antibody targeting asymmetric dimethylation at arginine 8 of Histone H3 (SWISS Q16695/P68431). Validated for WB, DB, IHC-P, and ELISA. Suitable for epigenetic and chromatin biology studies.

[Learn More](#)

Parameter	Value
Product Name	Asymmetric DiMethyl-Histone H3-R8 Rabbit mAb
Remarks/Alias	H3.4; H3/g; H3FT; H3t; HIST3H3; Histone H3; HIST1H3A; Asymmetric DiMethyl-Histone H3-R8
Species	Human
GeneID (human)	8290/8350
GeneID	8290/8350
Immunogen	A synthetic asymmetric dimethylated peptide around R8 of human Histone H3 (NP_003520.1)
Source	Rabbit
Category	Methylated Antibodies
Application	WB, DB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag