


Asymmetric Dimethyl-Histone H4-R3 (P62805) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa

Item Number: CM0000276



Introduction

Rabbit polyclonal antibody specific for asymmetric dimethylation at Arg3 of histone H4 (Swiss-Prot P62805). Validated for Western blotting, immunohistochemistry (paraffin), and ELISA. Predicted to react with human, mouse, rat, and other species, making it a valuable tool for epigenetic and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	Asymmetric DiMethyl-Histone H4-R3 Rabbit pAb
Remarks/Alias	H4; H4/n; H4C1; H4C2; H4C3; H4C4; H4C5; H4C6; H4C8; H4C9; H4F2; H4FN; FO108; H4-16; H4C11; H4C12; H4C13; H4C15; H4C16; HIST2H4; HIST2H4A; Asymmetric DiMethyl-Histone H4-R3
Species	Human
GeneID (Human)	H4
GeneID	8359
Immunogen	Synthetic peptide; A synthetic asymmetric dimethylated peptide around R3 of human HIST2H4A (NP_003529.1)
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P62805
Protein Weight	11kDa
Shipping	Ice bag