


Asymmetric Dimethyl-Histone H4-R3 Rabbit Pab | P62805 | Validated For Wb, If/Icc, Ihc-P, Elisa

Item Number: CM0000262



Introduction

Rabbit polyclonal antibody specific for asymmetric di-methylation of histone H4 at arginine 3 (P62805). Validated for WB, IF/ICC, IHC-P, and ELISA. Shows wide cross-reactivity with human, mouse, rat, and other species. Ideal for epigenetic and chromatin studies.

[Learn More](#)

Parameter	Value
Product Name	Asymmetric DiMethyl-Histone H4-R3 Rabbit pAb
Remarks/Alias	H4; H4/n; H4C1; H4C2; H4C3; H4C4; H4C5; H4C6; H4C8; H4C9; H4F2; H4FN; FO108; H4-16; H4C11; H4C12; H4C13; H4C15; H4C16; HIST2H4; HIST2H4A; Asymmetric DiMethyl-Histone H4-R3
Species	Human
GeneID (human)	H4
GeneID	8359
Immunogen	A synthetic asymmetric dimethylated peptide around R3 of human histone H4 (NP_003539.1)
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P62805
Protein Weight	11 kDa
Shipping	Ice bag