


Atad2 Rabbit Mab For Wb, Elisa - Q6PL18

Item Number: CM0003444



Introduction

ATAD2 Rabbit Monoclonal Antibody for Western Blot and ELISA, specific to human ATAD2 (Q6PL18). Ideal for research on breast cancer, estrogen receptor signaling, and transcriptional coactivation.

[Learn More](#)

Parameter	Value
Product Name	ATAD2 Rabbit mAb
Remarks/Alias	ANCCA; CT137; PRO2000; ATAD2
Species	Human
Gene ID (Human)	29028
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ATAD2 (Q6PL18)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q6PL18
Protein Weight	159kDa
Shipping	Ice bag