

Ate1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - 095260

Item Number: CM0007212



Introduction

Recombinant rabbit monoclonal antibody against human ATE1, validated for WB, IHC-P, and ELISA. Recognizes human, mouse, and rat ATE1. Ideal for studying N-terminal arginylation and ubiquitin-dependent degradation.

[Learn More](#)

Parameter	Value
Product Name	ATE1 Rabbit mAb
Species	Human
GeneID (Human)	11101
GeneID	11101
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human ATE1 (NP_001001976.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95260
Protein Weight	59kDa
Shipping	Ice bag