

Atf6 Rabbit Mab - P18850

Item Number: CM0007794

Introduction

Anti-ATF6 rabbit monoclonal antibody for WB, IP, IHC-P, and ELISA. Targets human ATF6, a key UPR transcription factor precursor embedded in the ER membrane. Reliable raw material for IVD development.

[Learn More](#)

Parameter	Value
Product Name	ATF6 Rabbit mAb
Remarks/Alias	ACHM7; ATF6A
Species	Human
GeneID	22926
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-194 of human ATF6 (NP_031374.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P18850
Protein Weight	75kDa
Shipping	Ice bag