

Atf6 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P18850

Item Number: CM0016689

Introduction

High-quality rabbit polyclonal antibody targeting ATF6 (ACHM7/ATF6A). Validated for WB, IHC-P, IF/ICC, and ELISA; cross-reacts with human, mouse, and rat. Ideal for ER stress and UPR research.

[Learn More](#)

Parameter	Specification
Product Name	ATF6 Rabbit pAb
Remarks/Alias	ACHM7; ATF6A; ATF6
Species	Human
GenelD (human)	22926
GenelD	22926
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18850
Protein Weight	75kDa
Shipping	Ice bag