

Atmin Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O43313

Item Number: CM0027087

Introduction

Rabbit polyclonal antibody targeting human ATMIN (ASCIZ/ZNF822). Validated for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Suitable for DNA damage and transcription research.

[Learn More](#)

Parameter	Value
Product Name	ATMIN Rabbit pAb
Remarks/Alias	ASCIZ; ZNF822; ATMIN
Species	Human
Gene ID (Human)	23300
Gene ID	23300
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 160-360 of human ATMIN (NP_056066.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43313
Protein Weight	88kDa
Shipping	Ice bag