

Atp1A2 Rabbit Pab For Wb, Elisa - P50993

Item Number: CM0014348

Introduction

High-quality rabbit polyclonal antibody targeting human ATP1A2 (P50993). Validated for WB and ELISA, cross-reacts with mouse and rat. Ideal for sodium pump research.

[Learn More](#)

Parameter	Value
Product Name	ATP1A2 Rabbit pAb
Remarks/Alias	FHM2; MHP2; DEE98; FARIMPD; ATP1A2
Species	Human
GeneID (Human)	477
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-70 of human ATP1A2 (NP_000693.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50993
Protein Weight	112kDa
Shipping	Ice bag