


Atp6V1A Rabbit Polyclonal Antibody For Wb, Elisa - P38606

Item Number: CM0013557



Introduction

Rabbit polyclonal antibody against human ATP6V1A, validated for Western blot and ELISA. Detects mouse protein. Targets the V-ATPase catalytic subunit for autophagy and lysosomal research.

[Learn More](#)

Parameter	Value
Product Name	ATP6V1A Rabbit pAb
Remarks/Alias	HO68; VA68; VPP2; Vma1; DEE93; ARCL2D; ATP6A1; IECEE3; ATP6V1A1; ATP6V1A
Species	Human
GeneID (Human)	523
GeneID	523
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 500-600 of human ATP6V1A (NP_001681.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P38606
Protein Weight	68kDa
Shipping	Ice bag