


Atr Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q13535

Item Number: CM0005019



Introduction

ATR rabbit monoclonal antibody validated for WB, IF/ICC, and ELISA in human samples. Targets the serine/threonine-protein kinase ATR, a critical DNA damage checkpoint regulator. UniProt Q13535.

[Learn More](#)

Parameter	Value
Product Name	ATR Rabbit mAb
Remarks/Alias	FRP1, MEC1, SCKL, FCTCS, SCKL1, ATR
Species	Human
GeneID	545
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 400-500 of human ATR (NP_001175.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q13535
Protein Weight	301kDa
Shipping	Ice bag