


Atr Rabbit Polyclonal Antibody For Wb, Elisa - Q13535

Item Number: CM0024774



Introduction

Rabbit polyclonal antibody targeting human ATR (serine/threonine-protein kinase ATR) for Western blot and ELISA applications. Validated for detection of ATR in DNA damage and checkpoint signaling research.

[Learn More](#)

Parameter	Value
Product Name	ATR Rabbit pAb
Remarks/Alias	FRP1; MEC1; SCKL; FCTCS; SCKL1; ATR
Species	Human
GeneID (human)	545
GeneID	545
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q13535
Protein Weight	301kDa
Shipping	Ice bag