

Aven Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nqs1

Item Number: CM0017679

Introduction

AVEN (PDCD12) rabbit polyclonal antibody for WB and ELISA. Detects human and rat AVEN, a cell death regulator that protects against Apaf-1-mediated apoptosis. Ideal for apoptosis studies in testis, ovary, thymus, and other tissues.

[Learn More](#)

Parameter	Value
Product Name	AVEN Rabbit pAb
Remarks/Alias	PDCD12; AVEN
Species	Human
Gene ID (Human)	57099
Gene ID	57099
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 203-362 of human AVEN (NP_065104.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9NQS1
Protein Weight	39kDa
Shipping	Ice bag