

Bag3 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - 095817

Item Number: CM0011209

Introduction

Rabbit monoclonal antibody against BAG3 (Bcl-2-associated athanogene 3), validated for WB, IHC-P, ELISA. Recognizes human and mouse BAG3; targets co-chaperone in autophagy, apoptosis, and stress response. Ideal for cancer and neurodegeneration studies.

[Learn More](#)

Parameter	Value
Product Name	BAG3 Rabbit mAb
Remarks/Alias	BIS; MFM6; BAG-3; CAIR-1; BAG3
Species	Human
GeneID	9531
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 50-150 of human BAG3 (095817).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	095817
Protein Weight	62kDa
Shipping	Ice bag