

Bank1 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Ndb2

Item Number: CM0015883

Introduction

BANK1 rabbit polyclonal antibody validated for WB and ELISA. Reacts with human, mouse, rat. Targets B-cell scaffold protein involved in BCR signaling and calcium mobilization. Suitable for immunology and autoimmunity research.

[Learn More](#)

Parameter	Value
Product Name	BANK1 Rabbit pAb
Remarks/Alias	BANK; BANK1
Species	Human
Gene ID	55024
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 516-785 of human BANK1 (NP_060405.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NDB2
Protein Weight	89kDa
Shipping	Ice bag