

Batf3 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Nr55

Item Number: CM0013622

Introduction

Rabbit polyclonal antibody against human BATF3, validated for Western blot and ELISA, with cross-reactivity to rat. Ideal for dendritic cell differentiation and immune response studies.

[Learn More](#)

Parameter	Value
Product Name	BATF3 Rabbit pAb
Remarks/Alias	JDP1; SNFT; JUNDM1; BATF3
Species	Human
GeneID (Human)	55509
GeneID	55509
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-127 of human BATF3 (NP_061134.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q9NR55
Protein Weight	14kDa
Shipping	Ice bag