


Baz1A Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nrl2

Item Number: CM0028502



Introduction

BAZ1A Rabbit Polyclonal Antibody validated for WB and ELISA. Reacts with human, mouse, rat. Targets the chromatin remodeling factor ACF1 involved in DNA replication and transcription. Ideal for research use.

[Learn More](#)

Parameter	Specification
Product Name	BAZ1A Rabbit pAb
Remarks/Alias	ACF1; WALp1; hACF1; WCRF180; BAZ1A
Species	Human
Gene ID (Human)	11177
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1347-1556 of human BAZ1A (NP_038476.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NRL2
Protein Weight	179kDa
Shipping	Ice bag