


Bcl10 Rabbit Polyclonal Antibody For Wb, Elisa - 095999

Item Number: CM0015202



Introduction

Rabbit polyclonal antibody targeting human Bcl10 (NP_003912.1, aa 1-118). Validated for WB and ELISA with cross-reactivity to human, mouse, and rat. Ideal for immune signaling and lymphoma research.

[Learn More](#)

Parameter	Specification
Product Name	Bcl10 Rabbit pAb
Remarks/Alias	CLAP; mE10; CIPER; IMD37; c-E10; CARMEN; Bcl10
Species	Human
GeneID (Human)	8915
GeneID	8915
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-118 of human Bcl10 (NP_003912.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95999
Protein Weight	26kDa
Shipping	Ice bag