

Bcl6 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P41182

Item Number: CM0029897



Introduction

BCL6 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat BCL6 protein (79kDa). BCL6 is a key transcriptional repressor in germinal center B-cells and lymphomas. Suitable for immunology and cancer research.

[Learn More](#)

Parameter	Value
Product Name	BCL6 Rabbit pAb
Remarks/Alias	BCL5; LAZ3; BCL6A; ZNF51; ZBTB27; BCL6
Species	Human
GeneID	604
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-280 of human BCL6 (NP_001124317.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P41182
Protein Weight	79kDa
Shipping	Ice bag