

Bcs1L Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y276

Item Number: CM0025144

Introduction

BCS1L Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, ELISA with cross-reactivity to human, mouse, rat. Targets mitochondrial chaperone BCS1L (48kDa) involved in complex III assembly; ideal for mitochondrial research and IVD development.

[Learn More](#)

Parameter	Value
Product Name	BCS1L Rabbit pAb
Remarks/Alias	BCS; BJS; PTD; BCS1; FLNMS; h-BCS; MC3DN1; h-BCS1; GRACILE; Hs.6719; BCS1L
Species	Human
GeneID (Human)	617
GeneID	617
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 33-240 of human BCS1L (NP_004319.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y276
Protein Weight	48kDa
Shipping	Ice bag