

Bgat Rabbit Monoclonal Antibody For Wb And Elisa - P16442

Item Number: CM0009225



Introduction

BGAT Rabbit Monoclonal Antibody validated for WB and ELISA, reacts with human, targeting the histo-blood group ABO system transferase (P16442). Ideal for blood group antigen research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	BGAT Rabbit mAb
Remarks/Alias	GTB; NAGAT; A3GALNT; A3GALT1; BGAT
Species	Human
GeneID (Human)	28
GeneID	28
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human BGAT (P16442).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P16442
Protein Weight	41kDa
Shipping	Ice bag