

Bgn Rabbit Polyclonal Antibody For Wb, Elisa - P21810

Item Number: CM0028907

Introduction

Rabbit polyclonal antibody targeting human Biglycan (BGN). Validated for WB and ELISA. Recombinant immunogen (aa 20-230). Ideal for extracellular matrix studies and collagen fiber assembly research.

[Learn More](#)

Parameter	Value
Product Name	BGN Rabbit pAb
Remarks/Alias	PGI; MRLS; DSPG1; PG-S1; SEMDX; SLRR1A; BGN
Species	Human
GeneID (human)	633
GeneID	633
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-230 of human BGN (NP_001702.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P21810
Protein Weight	42kDa
Shipping	Ice bag