


Bhmt Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q93088

Item Number: CM0007640



Introduction

BHMT Rabbit Monoclonal Antibody for WB, IHC-P, ELISA. Detects human and rat BHMT, 45kDa. Ideal for homocysteine metabolism and liver/kidney research.

[Learn More](#)

Parameter	Value
Product Name	BHMT Rabbit mAb
Remarks/Alias	BHMT1; HEL-S-61p; BHMT
Species	Human
Gene ID	635
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 306-406 of human BHMT. (NP_001704.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q93088
Protein Weight	45kDa
Shipping	Ice bag