

Bid Rabbit Pab - P55957

Item Number: CM0016224



Introduction

Bid Rabbit pAb (SWISS P55957) is a rabbit polyclonal antibody targeting human Bid. Suitable for WB, IF/ICC, ELISA, and cross-reactive with mouse and rat. Ideal for apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	Bid Rabbit pAb
Remarks/Alias	FP497; Bid
Species	Human
Gene ID	637
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-195 of human Bid (NP_001187.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55957
Protein Weight	22kDa
Shipping	Ice bag