

Biotin Anti-Cd4 Monoclonal Antibody For Fc - P01730

Item Number: CM0011911



Introduction

Biotin-conjugated rabbit monoclonal antibody against human and cynomolgus CD4, validated for flow cytometry. Ideal for T-helper cell identification, HIV receptor studies, and immune profiling.

[Learn More](#)

Parameter	Value
Product Name	Biotin Rabbit anti-Human/Monkey CD4 mAb
Remarks/Alias	T4; IMD79; Leu-3; OKT4D; CD4mut
Species	Human
Gene ID	920 (Human)
Immunogen	Recombinant protein of human CD4
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01730
Protein Weight	51kDa
Shipping	Ice bag