

Biotin Anti-Cd4 Monoclonal Antibody For Flow Cytometry - P01730

Item Number: CM0012617



Introduction

Biotin-conjugated rabbit monoclonal antibody targeting human and cynomolgus monkey CD4, validated for flow cytometry. Ideal for T-helper cell identification, immune response studies, and HIV research. Research-grade IVD raw material.

[Learn More](#)

Parameter	Value
Product Name	Biotin Rabbit anti-Human/Monkey CD4 mAb
Remarks/Alias	T4; IMD79; Leu-3; OKT4D; CD4mut
Species	Human
GeneID (Human)	920
GeneID	920
Immunogen	Recombinant protein of human CD4
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01730
Protein Weight	51 kDa
Shipping	Ice bag