


Biotin Anti-Cd8A Monoclonal Antibody For Flow Cytometry - P01732

Item Number: CM0011695



Introduction

Biotin-labeled rabbit monoclonal antibody targeting human and cynomolgus monkey CD8a, validated for flow cytometry. Ideal for T-cell immunophenotyping, immune checkpoint research, and cellular cytotoxicity assays.

[Learn More](#)

Parameter	Value
Product Name	Biotin Rabbit anti-Human/Monkey CD8a mAb
Remarks/Alias	CD8; p32; Leu2; CD8alpha
Species	Human
GeneID (Human)	925
GeneID	925
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-182 of human CD8A (NP_001759.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01732
Protein Weight	26 kDa
Shipping	Ice bag