

Biotin Anti-Hla-Dr Monoclonal Antibody For Flow Cytometry - P01903

Item Number: CM0012105



Introduction

Biotinylated rabbit monoclonal antibody targeting human/cynomolgus HLA-DR alpha chain (P01903), validated for flow cytometry. Ideal for MHC class II expression analysis on antigen-presenting cells and immune monitoring.

[Learn More](#)

Parameter	Value
Product Name	Biotin Rabbit anti-Human/Monkey HLA-DR mAb
Remarks/Alias	HLA-DRA1
Species	Human
GeneID (Human)	3122
GeneID	3122
Immunogen	Recombinant protein of human HLA-DR
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01903
Protein Weight	29kDa
Shipping	Ice bag