

Biotin Anti-Human Cd11B Monoclonal Antibody For Flow Cytometry - P11215

Item Number: CM0010890



Introduction

Biotin-conjugated mouse anti-human CD11b monoclonal antibody for flow cytometry. Recognizes human integrin alpha-M (MAC-1), expressed on monocytes and granulocytes. Suitable for cell surface staining and adhesion research.

[Learn More](#)

Parameter	Value
Product Name	Biotin Mouse anti-Human CD11b mAb
Remarks/Alias	CR3A, MO1A, CD11B, MAC-1, MAC1A, SLEB6
Species	Human
GeneID (Human)	3684
GeneID	3684
Immunogen	Synthetic peptide derived from human CD11b
Source	Mouse
Category	Monoclonal Antibodies
Application	FC (Flow Cytometry)
Cross Reactivity	Human
SWISS	P11215
Protein Weight	127 kDa
Shipping	Ice bag