

Biotin Anti-Mouse Cd117/C-Kit Monoclonal Antibody For Flow Cytometry - P05532

Item Number: CM0005953



Introduction

Biotin-conjugated rabbit monoclonal antibody targeting mouse CD117/c-Kit receptor tyrosine kinase. Validated for flow cytometry, detects extracellular domain. Ideal for stem cell biology, hematopoiesis, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Biotin Rabbit anti-Mouse CD117/c-Kit mAb
Remarks/Alias	W; Bs; Fdc; Ssm; SCO1; SCO5; SOW3; CD117; c-KIT; Tr-kit; Gsfsc01; Gsfsc05; Gsfscow3
Species	Mouse
GeneID (Mouse)	16590
GeneID (Human)	3815
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-523 of mouse CD117/c-Kit (NP_001116205.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P05532
Protein Weight	22kDa/108kDa/109kDa
Shipping	Ice bag