

Biotin Anti-Mouse Cd24 Monoclonal Antibody For Flow Cytometry - P24807

Item Number: CM0005849



Introduction

Biotin rabbit anti-mouse CD24 monoclonal antibody for flow cytometry. Targets a GPI-anchored glycoprotein involved in thymocyte development, B-cell modulation, and immune suppression. Validated in mouse.

[Learn More](#)

Parameter	Value
Product Name	Biotin Rabbit anti-Mouse CD24 mAb
Remarks/Alias	HSA; Cd24; Ly-52; nectadrin
Species	Mouse
Gene ID (Human)	100133941
Gene ID	12484
Immunogen	Recombinant protein (amino acids 27-53 of mouse CD24, NP_033976.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P24807
Protein Weight	8kDa
Shipping	Ice bag