

Biotin Anti-Mouse Tcr Γ/Δ Monoclonal Antibody - Tcr Γ/Δ

Item Number: CM0011963

Introduction

Biotin-conjugated Armenian Hamster monoclonal antibody targeting mouse TCR γ/δ , optimized for flow cytometry. Ideal for studying gamma delta T cell subsets in immunological research.

[Learn More](#)

Parameter	Value
Product Name	Biotin anti-Mouse TCR γ/δ mAb
Remarks/Alias	Tcrdelta; TCRGV1S1; TCRGV2S1; TCRGV3S1; TCRGV5S3
Species	Mouse
Gene ID	110066/110067
Immunogen	Recombinant protein of mouse TCR γ/δ
Source	Armenian Hamster
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
Shipping	Ice bag