

Biotin Rabbit Anti-Human Cd34 Monoclonal Antibody For Flow Cytometry - P28906

Item Number: CM0011191

Introduction

Biotinylated rabbit monoclonal antibody targeting human CD34 (P28906) for flow cytometry. Validated on human hematopoietic progenitor cells. Reliable for stem cell research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	Biotin Rabbit anti-Human CD34 mAb
Remarks/Alias	CD34;CD34 molecule;GIG3
Species	Human
GenelD (Human)	947
GenelD	947
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 170-290 of human CD34 (NP_001020280.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P28906
Protein Weight	35 kDa/41 kDa
Shipping	Ice bag