

Biotin Rabbit Anti-Human Cd34 Monoclonal Antibody For Flow Cytometry - P28906

Item Number: CM0005083

Introduction

Biotin-conjugated rabbit monoclonal antibody targeting human CD34, a hematopoietic progenitor cell marker. Validated for flow cytometry. Cross-reacts with human. Suitable for stem cell research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	Biotin Rabbit anti-Human CD34 mAb
Remarks/Alias	CD34; CD34 molecule; GIG3
Species	Human
Gene ID (Human)	947
Gene ID	947
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 170-290 of human CD34 (NP_001020280.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P28906
Protein Weight	35kDa/41kDa
Shipping	Ice bag